

CrushedCRMs™ for PAL (Testing Recoverable Au)

The use of pulverised CRMs in PAL can be problematic.

Our CrushedCRM™ - Pre-Mixed product restores confidence in the PAL technique by unlocking the following solutions:

- CrushedCRM™ - Pre-Mixed offers a more precise solution based on CRM mass and grade.
- CrushedCRM™ - Pre-Mixed are packaged at the aliquot size. They are a single-shot, ready-to-use product.
- CrushedCRM™ - Pre-Mixed test the pulverisation step in the process, and the gold is leached over time, unlike pulverised gold CRMs, which will leach very quickly as the gold is already liberated.
- CrushedCRM™ - Pre-Mixed can be submitted blindly to the lab from the field, unlike pulverised gold CRMs.
- CrushedCRM™ - Pre-Mixed are designed and certified for the PAL method. Some gold will be retained in the tails, allowing for a more complete quality control program to be established from a single source.
- CrushedCRM™ - Pre-Mixed are also certified by PhotonAssay™ and Fire Assay for total gold content.

The key benefits of using CrushedCRM™ for PAL are:

- CrushedCRMs™ are pulverised and leached in the method like a routine sample.
- CrushedCRMs™ are certified for PAL at the method aliquot mass.