

Crushed CRM™ GOLD

PRE-MIXED

CrushedCRM™ – Pre-Mixed is a ready-to-use single-shot gold CRM certified for PAL and PhotonAssay™ at the method aliquot mass.

Product Benefits

- Certification for the aliquot mass method.
- Tests the entire sample preparation process.
- Leaches over time like geological samples.
- Requires pulverising during the test.
- QC confidence over time.
- Sample-to-sample variance similar to Pulverised Gold.
- Easy to swap to CrushedCRM™ from Pulverised Gold.
- Long-term, guaranteed product supply.
- PAL & PhotonAssay™ certifications provide trust and confidence in your results.

Product Features

- Crushed ~3mm rock chip.
- Certified by PAL and PhotonAssay™ method through international round robin, treated the same way as the analytical method.
- Within unit heterogeneity to identify sample preparation issues.
- Between unit homogeneity similar to pulverised CRM.
- Large production batch sizes (multi-ton).

Product Usage

CrushedCRM™ - Pre-Mixed CRMs are suitable for use:

- PhotonAssay™.
- PAL.
- Fire Assay.

They are available in:

- 500g units. Ready to use. Suitable for PhotonAssay™ and PAL. The certified 'value' and 'uncertainty' are certified for the PAL and aliquot mass methods.
- 2kg Pack. This tests the uncertainty of sample preparation (crushed, split, pulverised). It is for use in the Fire Assay.

Product Grades

It is available in the following grades from 0.1 – 3.0.